



## Product description

Order number	70.3030.100
Product description	Pipette tip, working volume: 200 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack

## Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	200 µl
Filter	no
Coding colour	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 200 µl, transparent, PCR Performance**  
**Tested, 480 piece(s)/StackPack**

**Page 2**

**Size**

Max. volume	200 µl
Length of product	51 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	1920
Type of smallest subpackaging	StackPack
Piece(s) / inner box	1920
Piece (s) / outer case	1920
Piece(s) / pallet	291840
Depth of case	290 mm
Width of case	157 mm
Height of case	225 mm
Case volume	0.0102 cbm
Weight of product	0.0006 kg
Weight of case	1.69 kg
EAN of inner box	4038917845960
EAN case	4038917845977

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-23

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**