




Product Specification
Filter tip, 20 µl, transparent, Biosphere® plus,
Low retention, 96 piece(s)/SingleRefill

Page 1



Product description

Order number	70.3030.385
Product description	Filter tip, working volume: 20 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/SingleRefill

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	20 µl
Filter	yes
Coding colour	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Filter tip, 20 µl, transparent, Biosphere® plus,
Low retention, 96 piece(s)/SingleRefill**

Page 2

Size

Max. volume	20 µl
Length of product	51 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	SingleRefill
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	115200
Depth of box	291 mm
Width of box	105 mm
Height of box	148 mm
Depth of case	437 mm
Width of case	307 mm
Height of case	170 mm
Case volume	0.0228 cbm
Weight of product	0.0008 kg
Weight of case	5.16 kg
EAN of inner box	4038917847070
EAN case	4038917847087

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT