



**Product description**

|                     |                                                                                                                                                                                                                                                            |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 70.3031.110                                                                                                                                                                                                                                                |
| Product description | Pipette tip, working volume: 250 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack |

**Product characteristics**

|                |                                                                                                         |
|----------------|---------------------------------------------------------------------------------------------------------|
| Graduation     | yes                                                                                                     |
| Suitable for   | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design |
| Volume of work | 250 µl                                                                                                  |
| Filter         | no                                                                                                      |
| Coding colour  | yellow                                                                                                  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Pipette tip, 250 µl, transparent, PCR Performance Tested, Low retention, 480 piece(s)/StackPack**

**Page 2**

### Size

|                   |        |
|-------------------|--------|
| Max. volume       | 250 µl |
| Length of product | 51 mm  |

### Material & colours

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

### Purity & certification

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Purity standard  | PCR Performance Tested        |
| Batched          | yes                           |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1920          |
| Type of smallest subpackaging | StackPack     |
| Piece(s) / inner box          | 1920          |
| Piece (s) / outer case        | 1920          |
| Piece(s) / pallet             | 291840        |
| Depth of case                 | 290 mm        |
| Width of case                 | 157 mm        |
| Height of case                | 225 mm        |
| Case volume                   | 0.0102 cbm    |
| Weight of product             | 0.0003 kg     |
| Weight of case                | 1.69 kg       |
| EAN of inner box              | 4038917845533 |
| EAN case                      | 4038917845540 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**