



Product description

| | |
|---------------------|--|
| Order number | 70.3040.305 |
| Product description | Pipette tip, working volume: 300 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/SingleRefill |

Product characteristics

| | |
|------------------------|---|
| Graduation | yes |
| Suitable for | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design |
| Volume of work | 300 µl |
| Filter | no |
| Coding colour | orange |
| Autoclavable (max. °C) | 121 °C |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Pipette tip, 300 µl, transparent, Biosphere® plus,
96 piece(s)/SingleRefill**

Page 2

Size

| | |
|-------------------|--------|
| Max. volume | 300 µl |
| Diameter | 7 mm |
| Length of product | 60 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|---------------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| Certificate of conformity | yes |
| CE certified | CE |
| Purity standard | Biosphere® plus |
| Sterilisation | Gas sterilisation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1920 |
| Type of smallest subpackaging | SingleRefill |
| Piece(s) / inner box | 480 |
| Piece (s) / outer case | 1920 |
| Piece(s) / pallet | 103680 |
| Depth of box | 333 mm |
| Width of box | 105 mm |
| Height of box | 148 mm |
| Depth of case | 437 mm |
| Width of case | 351 mm |
| Height of case | 172 mm |
| Case volume | 0.0264 cbm |
| Weight of product | 0.0009 kg |
| Weight of case | 2.74 kg |
| EAN of inner box | 4038917433136 |
| EAN case | 4038917433129 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

