



Product description

Order number	70.3050.030
Product description	Pipette tip, working volume: 1,000 µl, black, PCR Performance Tested, 500 piece(s)/bag

Product characteristics

Graduation	no
Volume of work	1.000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipette tip, 1,000 µl, black, PCR Performance
Tested, 500 piece(s)/bag

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	black

Purity & certification

Product category	no medical device IVD
Certificate of conformity	yes
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0009 kg
Weight of case	5.24 kg
EAN of inner box	4038917434614
EAN case	4038917434607

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT