




Product description

Order number	70.3050.100
Product description	Pipette tip, working volume: 1,000 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	1.000 µl
Filter	no
Coding colour	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipette tip, 1,000 µl, transparent, PCR
Performance Tested, 480 piece(s)/StackPack

Page 2

Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	StackPack
Piece(s) / inner box	1920
Piece (s) / outer case	1920
Piece(s) / pallet	149760
Depth of case	303 mm
Width of case	295 mm
Height of case	225 mm
Case volume	0.0201 cbm
Weight of product	0.0012 kg
Weight of case	2.82 kg
EAN of inner box	4038917424844
EAN case	4038917424837

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT