




Product Specification
**Pipette tip, 1,250 µl, transparent, Biosphere®
plus, 96 piece(s)/SingleRefill**

Page 1



Product description

Order number	70.3060.305
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/SingleRefill

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design
Volume of work	1.250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipette tip, 1,250 µl, transparent, Biosphere®
plus, 96 piece(s)/SingleRefill

Page 2

Size

Max. volume	1250 µl
Diameter	9 mm
Length of product	90 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	SingleRefill
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	69120
Depth of box	495 mm
Width of box	105 mm
Height of box	148 mm
Depth of case	512 mm
Width of case	325 mm
Height of case	237 mm
Case volume	0.0394 cbm
Weight of product	0.0016 kg
Weight of case	4.38 kg
EAN of inner box	4038917432474
EAN case	4038917432467

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT