



### Product description

Order number	90.3100.010
Product description	Sarpette® M, single-channel, variable, 0.5 µl - 10 µl, grey, 1 piece(s)/case

### Product characteristics

Volume of work	0,5 µl - 10 µl
Coding colour	grey
Autoclavable (max. °C)	121 °C

### Material & colours

Colour of product	grey
-------------------	------

### Purity & certification

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	315 mm
Width of box	90 mm
Height of box	30 mm
Weight of product	0.088 kg
EAN of inner box	4038917936729
EAN case	4038917936736

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

