



Product Specification
**Pipetting aids, Rotatable, storage capacity: 6, for
Single-channel pipettes up to 1,000 µl, 1
piece(s)/case**

Page 1



Product description

Order number	95.1760.200
Product description	Pipetting aids, Rotatable, storage capacity: 6, for Single-channel pipettes up to 1,000 µl, 1 piece(s)/case

Product characteristics

Suitable for	Single-channel pipettes up to 1,000 µl
Storage capacity	6

Material & colours

Colour of product	white/red
-------------------	-----------

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipetting aids, Rotatable, storage capacity: 6, for
Single-channel pipettes up to 1,000 µl, 1
piece(s)/case

Page 2

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	4
Piece(s) / pallet	144
Depth of box	230 mm
Width of box	150 mm
Height of box	150 mm
Depth of case	360 mm
Width of case	330 mm
Height of case	250 mm
Case volume	0.0297 cbm
Weight of product	0.68 kg
Weight of case	3.16 kg
EAN of inner box	4038917937863
EAN case	4038917937870

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT