




## Product Specification

**Serological pipette, with tip, plugged, 1 ml, sterile,  
non-pyrogenic/endotoxin-free, non-cytotoxic, 1  
piece(s)/blister**

**Page 1**



### Product description

Order number	86.1251.001
Product description	Serological pipette, with tip, plugged, 1 ml, coding colour: yellow, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

### Product characteristics

Colour of print/label	yellow
Graduation	yes
Design	with tip, plugged
Volume of work	1 ml
Coding colour	yellow

### Size

Diameter	4 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Serological pipette, with tip, plugged, 1 ml, sterile,  
non-pyrogenic/endotoxin-free, non-cytotoxic, 1  
piece(s)/blister**

**Page 2**

### Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

### Purity & certification

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	54000
Depth of case	420 mm
Width of case	250 mm
Height of case	290 mm
Case volume	0.0323 cbm
Weight of product	0.0044 kg
Weight of case	5.01 kg
EAN of inner box	4038917127158
EAN case	4038917021067

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**