




Product Specification
Serological pipette, with tip, plugged, 1 ml, sterile,
25 piece(s)/bag

Page 1



Product description

Order number	86.1251.025
Product description	Serological pipette, with tip, plugged, 1 ml, coding colour: yellow, material: PS, sterile, 25 piece(s)/bag

Product characteristics

Colour of print/label	yellow
Graduation	yes
Design	with tip, plugged
Volume of work	1 ml
Coding colour	yellow

Size

Diameter	4 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Serological pipette, with tip, plugged, 1 ml, sterile,
25 piece(s)/bag

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	1000
Piece(s) / pallet	108000
Depth of case	413 mm
Width of case	250 mm
Height of case	150 mm
Case volume	0.0155 cbm
Weight of product	0.0033 kg
Weight of case	3.94 kg
EAN of inner box	4038917127165
EAN case	4038917021074

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-24

This document was prepared by EDP support and is valid without signature.



SARSTEDT