




Product Specification
**Serological pipette, without tip, plugged, 10 ml,
sterile, 10 piece(s)/bag**

Page 1



Product description

Order number	86.1688.010
Product description	Serological pipette, without tip, plugged, 10 ml, material: PS, sterile, 10 piece(s)/bag

Product characteristics

Graduation	yes
Design	without tip, plugged
Volume of work	10 ml

Size

Diameter	10 mm
Length of product	335 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Serological pipette, without tip, plugged, 10 ml,
sterile, 10 piece(s)/bag**

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	500
Piece(s) / pallet	27000
Depth of case	413 mm
Width of case	250 mm
Height of case	270 mm
Case volume	0.0279 cbm
Weight of product	0.0085 kg
Weight of case	4.97 kg
EAN of inner box	4038917133081
EAN case	4038917024693

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

