




Product Specification
Transfer pipette, 1 ml, (LxW): 115 x 10 mm, LD-PE,
transparent

Page 1



Product description

Order number	86.1170.300
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 115 x 10 mm, material: LD-PE, transparent, droplet size: 40 µl - 50 µl, 500 piece(s)/case

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	1 ml
Width of product	10 mm
Length of product	115 mm
Drop size	40 µl - 50 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Transfer pipette, 1 ml, (LxW): 115 x 10 mm, LD-PE,
transparent

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	4000
Piece(s) / pallet	156000
Case volume	0.0401 cbm
Weight of product	0.001 kg
Weight of case	4.55 kg
EAN of inner box	4038917378352
EAN case	4038917378345

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

