




**Product Specification**  
**Transfer pipette, 1 ml, (LxW): 115 x 10 mm, LD-PE,**  
**transparent**

**Page 1**



**Product description**

Order number	86.1170
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 115 x 10 mm, material: LD-PE, transparent, droplet size: 40 µl - 50 µl, injected graduation up to 0.3 ml, 500 piece(s)/case

**Product characteristics**

Graduation	yes
------------	-----

**Size**

Sample volume	1 ml
Width of product	10 mm
Length of product	115 mm
Drop size	40 µl - 50 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 1 ml, (LxW): 115 x 10 mm, LD-PE,  
transparent**

**Page 2**

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	168000
Depth of box	261 mm
Width of box	126 mm
Height of box	121 mm
Depth of case	533 mm
Width of case	263 mm
Height of case	140 mm
Case volume	0.0196 cbm
Weight of product	0.0007 kg
Weight of case	1.89 kg
EAN of inner box	4038917085564
EAN case	4038917004206

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-09

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**