




**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,**  
**LD-PE, transparent**

**Page 1**



**Product description**

Order number	86.1171.001
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 15 mm, material: LD-PE, transparent, droplet size: 35 µl - 55 µl, injected graduation up to 3 ml, sterile, 1 piece(s)/blister

**Product characteristics**

Graduation	yes
------------	-----

**Size**

Sample volume	3.5 ml
Width of product	15 mm
Length of product	155 mm
Drop size	35 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,**  
**LD-PE, transparent**

**Page 2**

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	840
Type of smallest subpackaging	blister
Piece(s) / inner box	42
Piece (s) / outer case	840
Piece(s) / pallet	30240
Depth of case	578 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0415 cbm
Weight of product	0.0015 kg
Weight of case	2.6 kg
EAN of inner box	4038917126267
EAN case	4038917020541

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**