




Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1171.020
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 15 mm, material: LD-PE, transparent, droplet size: 35 µl - 55 µl, sterile, 20 piece(s)/bag

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3.5 ml
Width of product	15 mm
Length of product	155 mm
Drop size	35 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,
LD-PE, transparent

Page 2

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	4000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	4000
Piece(s) / pallet	64000
Case volume	0.1359 cbm
Weight of product	0.002 kg
Weight of case	10.18 kg
EAN of inner box	4038917318280
EAN case	4038917318273

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT