




**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,**  
**LD-PE, transparent**

**Page 1**



**Product description**

Order number	86.1171.300
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 15 mm, material: LD-PE, transparent, droplet size: 35 µl - 55 µl, 500 piece(s)/case

**Product characteristics**

Graduation	yes
------------	-----

**Size**

Sample volume	3.5 ml
Width of product	15 mm
Length of product	155 mm
Drop size	35 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,**  
**LD-PE, transparent**

**Page 2**

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	56000
Case volume	0.0994 cbm
Weight of product	0.0017 kg
Weight of case	9.41 kg
EAN of inner box	4038917378383
EAN case	4038917378376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**