



Product Specification
Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1172.020
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, sterile, 20 piece(s)/bag

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	156 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,
LD-PE, transparent

Page 2

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 4000
Type of smallest bag
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 4000
Piece(s) / pallet 64000
Case volume 0.1309 cbm
Weight of product 0.0017 kg
Weight of case 8.64 kg
EAN of inner box 4038917318310
EAN case 4038917318303

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT