



Product Specification
Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1172.300
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	156 mm
Drop size	30 µl - 45 µl

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	56000
Case volume	0.0994 cbm
Weight of product	0.0018 kg
Weight of case	10 kg
EAN of inner box	4038917378413
EAN case	4038917378406

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

