




**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 12.5 mm,**  
**LD-PE, transparent**

**Page 1**



**Product description**

Order number	86.1173.001
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, sterile, 1 piece(s)/blister

**Product characteristics**

Graduation	no
------------	----

**Size**

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	155 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	4000
Type of smallest subpackaging	blister
Piece(s) / inner box	500
Piece (s) / outer case	4000
Piece(s) / pallet	64000
Case volume	0.1309 cbm
Weight of product	0.0023 kg
Weight of case	11.23 kg
EAN of inner box	4038917318341
EAN case	4038917318334

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

