




**Product Specification**  
**Transfer pipette, 6 ml, (LxW): 152 x 15 mm, LD-PE,**  
**transparent**

**Page 1**



**Product description**

Order number	86.1174
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 152 x 15 mm, material: LD-PE, transparent, droplet size: 40 µl - 55 µl, 400 piece(s)/case

**Product characteristics**

Graduation	no
------------	----

**Size**

Sample volume	6 ml
Width of product	15 mm
Length of product	152 mm
Drop size	40 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 6 ml, (LxW): 152 x 15 mm, LD-PE,**  
**transparent**

**Page 2**

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	1600
Type of smallest subpackaging	case
Piece(s) / inner box	400
Piece (s) / outer case	1600
Piece(s) / pallet	44800
Depth of box	261 mm
Width of box	161 mm
Height of box	231 mm
Depth of case	533 mm
Width of case	333 mm
Height of case	250 mm
Case volume	0.0444 cbm
Weight of product	0.0017 kg
Weight of case	3.58 kg
EAN of inner box	4038917085601
EAN case	4038917004244

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-06

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**