




**Product Specification**  
**Transfer pipette, 6 ml, (LxW): 146 x 15 mm, LD-PE,**  
**transparent**

**Page 1**



**Product description**

Order number	86.1175.300
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 146 x 15 mm, material: LD-PE, transparent, droplet size: 20 µl - 30 µl

**Product characteristics**

Graduation	no
------------	----

**Size**

Sample volume	6 ml
Width of product	15 mm
Length of product	146 mm
Drop size	20 µl - 30 µl

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	400
Piece(s) / inner box	400
Piece (s) / outer case	4000
Piece(s) / pallet	44800
Case volume	0.0994 cbm
Weight of product	0.0018 kg
Weight of case	7.95 kg
EAN of inner box	4038917378536
EAN case	4038917378529

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

