



Product Specification
Transfer pipette, 1 ml, (LxW): 87 x 10 mm, LD-PE,
transparent

Page 1



Product description

Order number	86.1179.300
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 87 x 10 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl

Product characteristics

Graduation	no
------------	----

Size

Sample volume	1 ml
Width of product	10 mm
Length of product	87 mm
Drop size	30 µl - 45 µl

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Transfer pipette, 1 ml, (LxW): 87 x 10 mm, LD-PE,
transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 500
Piece(s) / inner box 500
Piece (s) / outer case 4000
Piece(s) / pallet 234000
Case volume 0.041 cbm
Weight of product 0.001 kg
Weight of case 4.55 kg
EAN of inner box 4038917378598
EAN case 4038917378581

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

