



**Product Specification**  
**DishRack, height: 370 mm, red, for 88 petri dishes with 92 mm Ø**

**Page 1**



**Product description**

Order number	93.1646.001
Product description	DishRack, height: 370 mm, material: PC, red, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case

**Product characteristics**

Storage capacity	88
------------------	----

**Size**

Height of product	370 mm
-------------------	--------

**Material & colours**

Product material	Polycarbonate (PC)
Colour of product	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category                      no medical device | IVD  
Batched                                    no

**Packaging**

Minimum order qty.                    1  
Type of smallest                        case  
subpackaging  
Piece(s) / inner box                    1  
Piece (s) / outer case                  1  
Piece(s) / pallet                        108  
Depth of case                            383 mm  
Width of case                            193 mm  
Height of case                            190 mm  
Case volume                              0.014 cbm  
Weight of product                       0.536 kg  
Weight of case                           0.8 kg  
EAN of inner box                        4038917132558  
EAN case                                  4038917030465

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

