



**Product Specification**  
**DishRack, height: 370 mm, blue, for 88 petri dishes with 92 mm Ø**

**Page 1**



**Product description**

Order number 93.1646.003  
Product description DishRack, height: 370 mm, material: PC, blue, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case

**Product characteristics**

Storage capacity 88

**Size**

Height of product 370 mm

**Material & colours**

Product material Polycarbonate (PC)  
Colour of product blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**DishRack, height: 370 mm, blue, for 88 petri dishes with 92 mm Ø**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      1  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      1  
Piece (s) / outer case                      1  
Piece(s) / pallet                              108  
Depth of case                              383 mm  
Width of case                              193 mm  
Height of case                              190 mm  
Case volume                              0.014 cbm  
Weight of product                              0.536 kg  
Weight of case                              0.8 kg  
EAN of inner box                              4038917132572  
EAN case                                      4038917030489

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

