



**Product Specification**  
**DishRack, height: 370 mm, green, for 88 petri  
dishes with 92 mm Ø**

**Page 1**



**Product description**

Order number	93.1646.004
Product description	DishRack, height: 370 mm, material: PC, green, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case

**Product characteristics**

Storage capacity	88
------------------	----

**Size**

Height of product	370 mm
-------------------	--------

**Material & colours**

Product material	Polycarbonate (PC)
Colour of product	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**DishRack, height: 370 mm, green, for 88 petri  
dishes with 92 mm Ø**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      1  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      1  
Piece (s) / outer case                      1  
Piece(s) / pallet                              108  
Depth of case                              383 mm  
Width of case                                193 mm  
Height of case                               190 mm  
Case volume                                0.014 cbm  
Weight of product                        1.086 kg  
Weight of case                              1.35 kg  
EAN of inner box                        4038917132589  
EAN case                                    4038917030496

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

