



## Product Specification

DishRack, height: 240 mm, blue, for 52 petri dishes with 92 mm Ø

Page 1



### Product description

Order number 93.1647.003  
Product description DishRack, height: 240 mm, material: PC, blue, for 52 petri dishes with 92 mm Ø, 1 piece(s)/case

### Product characteristics

Storage capacity 52

#### Size

Height of product 240 mm

### Material & colours

Product material Polycarbonate (PC)  
Colour of product blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**DishRack, height: 240 mm, blue, for 52 petri dishes with 92 mm Ø**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      1  
Type of smallest subpackaging                      case  
Piece(s) / inner box                      1  
Piece (s) / outer case                      1  
Piece(s) / pallet                      162  
Depth of case                      253 mm  
Width of case                      193 mm  
Height of case                      190 mm  
Case volume                      0.0093 cbm  
Weight of product                      0.617 kg  
Weight of case                      0.8 kg  
EAN of inner box                      4038917132626  
EAN case                      4038917030533

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

