



Product Specification
Semi-micro cuvette, 3 ml, (HxW): 45 x 12 mm, PS,
transparent, optical sides: 2

Page 1



Product description

Order number	67.746
Product description	Semi-micro cuvette, max. working volume: 3 ml, material: PS, transparent, (HxW): 45 x 12 mm, optical sides: 2, path length: 10 mm, 500 piece(s)/bag

Product characteristics

Volume of work	1 ml - 3 ml
Number of optical sides	2
Path length	10 mm
Wavelength in the cuvette area	330 nm

Size

Max. volume	3 ml
Width of product	12 mm
Height of product	45 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Semi-micro cuvette, 3 ml, (HxW): 45 x 12 mm, PS,
transparent, optical sides: 2

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	96000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0023 kg
Weight of case	5.02 kg
EAN of inner box	4038917083928
EAN case	4038917002806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT