




Product Specification
**Petri dish, 92 x 16 mm, transparent, with
ventilation cams**

Page 1



Product description

Order number	82.1473.001
Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

Size

Max. volume	80 ml
Diameter	92 mm
Height of product	16 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Petri dish, 92 x 16 mm, transparent, with
ventilation cams**

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Purity standard	sterile
Sterilisation	Gamma irradiation
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	480
Type of smallest subpackaging	tube bag
Piece(s) / inner box	480
Piece (s) / outer case	480
Piece(s) / pallet	9600
Depth of case	584 mm
Width of case	387 mm
Height of case	360 mm
Case volume	0.0814 cbm
Weight of product	0.016 kg
Weight of case	8.79 kg
EAN of inner box	4038917129329
EAN case	4038917022125

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT