



Product description

Order number	82.1473.020
Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, red, with ventilation cams, 20 piece(s)/tube bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

Size

Max. volume	80 ml
Diameter	92 mm
Height of product	16 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Latex-free	yes
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	480
Type of smallest subpackaging	tube bag
Piece(s) / inner box	480
Piece (s) / outer case	480
Piece(s) / pallet	9600
Depth of case	584 mm
Width of case	387 mm
Height of case	360 mm
Case volume	0.0814 cbm
Weight of product	0.0169 kg
Weight of case	9.25 kg
EAN of inner box	4038917129343
EAN case	4038917022149

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

