



**Product Specification**  
**Petri dish, 92 x 16 mm, green, with ventilation**  
**cams**

**Page 1**



**Product description**

Order number	82.1473.080
Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, green, with ventilation cams, 20 piece(s)/tube bag

**Product characteristics**

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

**Size**

Max. volume	80 ml
Diameter	92 mm
Height of product	16 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polystyrene (PS)
Colour of product	green

**Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Safety data sheet	no
Batched	yes

**Packaging**

Minimum order qty.	480
Type of smallest subpackaging	tube bag
Piece(s) / inner box	480
Piece (s) / outer case	480
Piece(s) / pallet	9600
Depth of case	584 mm
Width of case	387 mm
Height of case	360 mm
Case volume	0.0814 cbm
Weight of product	0.0169 kg
Weight of case	9.25 kg
EAN of inner box	4038917129404
EAN case	4038917022170

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

