



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 82.9923.426  |
| Product description | Petri dish, (ØxH): 54.65 x 14.7 mm, material: PS, transparent, without ventilation cams, 10 piece(s)/bag |

### **Product characteristics**

|                  |    |
|------------------|----|
| Ventilation cams | no |
|------------------|----|

### **Size**

|                   |          |
|-------------------|----------|
| Diameter          | 54.65 mm |
| Height of product | 14.7 mm  |

### **Material & colours**

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | transparent      |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

|                   |     |
|-------------------|-----|
| Latex-free        | yes |
| Pyrogen-free      | no  |
| Safety data sheet | no  |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 10            |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 40000         |
| Depth of case                 | 604 mm        |
| Width of case                 | 264 mm        |
| Height of case                | 373 mm        |
| Case volume                   | 0.06 cbm      |
| Weight of case                | 8.2 kg        |
| EAN of inner box              | 4038917154208 |
| EAN case                      | 4038917039772 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

