



Product description

Order number	82.9923.427
Product description	Petri dish, (ØxH): 54.65 x 14.7 mm, material: PS, transparent, without ventilation cams, sterile, 10 piece(s)/bag

Product characteristics

Ventilation cams	no
------------------	----

Size

Diameter	54.65 mm
Height of product	14.7 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	1000
Piece(s) / pallet	40000
Depth of case	604 mm
Width of case	264 mm
Height of case	373 mm
Case volume	0.06 cbm
Weight of case	8.2 kg
EAN of inner box	4038917154215
EAN case	4038917039789

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

