



Product description

| | |
|---------------------|---|
| Order number | 80.1301 |
| Product description | Swabs, in the tube, no medium, (LxØ): 175 x 12 mm, stem material: PS, stem length: 131 mm, swab material: viscose, sterile, 100 piece(s)/case |

Product characteristics

| | |
|----------------|-------------|
| Closure type | push cap |
| Design | in the tube |
| Type of medium | no medium |

Size

| | |
|----------------------------------|--------|
| Diameter | 12 mm |
| Length of product | 175 mm |
| Length including cap | 175 mm |
| Length excluding cap | 155 mm |
| Length including cap and plunger | 166 mm |
| Length of stem | 131 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | white |
| Swab material | viscose |
| Stem material | Polystyrene (PS) |

Purity & certification

| | |
|-------------------|-------------------|
| Product category | Medical device |
| CE certified | CE supplier |
| Purity standard | sterile |
| Sterilisation | Gamma irradiation |
| Safety data sheet | no |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 32000 |
| Depth of case | 590 mm |
| Width of case | 220 mm |
| Height of case | 300 mm |
| Case volume | 0.0389 cbm |
| Weight of product | 0.0047 kg |
| Weight of case | 5.94 kg |
| EAN of inner box | 4038917086103 |
| EAN case | 4038917234184 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

