




## Product Specification

Block rotor, for 40 tubes up to 8.5 mm Ø  
(S-Sedivette®), for SARMIX® M 2000

Page 1



### Product description

Order number	92.180.615
Product description	Block rotor, for 40 tubes up to 8.5 mm Ø (S-Sedivette®), for SARMIX® M 2000, 1 piece(s)/case

### Product characteristics

Suitable for	SARMIX® M 2000
Storage capacity	40

### Purity & certification

Product category	no medical device   IVD
Batched	no

### Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Block rotor, for 40 tubes up to 8.5 mm Ø  
(S-Sedivette®), for SARMIX® M 2000**

**Page 2**

Case volume	0.0122 cbm
Weight of product	0.2509 kg
Weight of case	2.71 kg
EAN of inner box	4038917102032
EAN case	4038917074315

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-20

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**