



Product description

Order number	92.180.618
Product description	Block rotor, for 36 tubes up to 15 mm diameter, for SARMIX® M 2000, 1 piece(s)/case

Product characteristics

Suitable for	SARMIX® M 2000
Storage capacity	36

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Block rotor, for 36 tubes up to 15 mm diameter,
for SARMIX® M 2000

Page 2

Depth of case	350 mm
Width of case	240 mm
Height of case	150 mm
Case volume	0.0126 cbm
Weight of product	0.2108 kg
Weight of case	2.31 kg
EAN of inner box	4038917309561
EAN case	4038917309554

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT