



### **Product description**

|                     |                                |
|---------------------|--------------------------------|
| Order number        | 92.181.611                     |
| Product description | Rotor, for 22 collection tubes |

### **Product characteristics**

|                  |    |
|------------------|----|
| Storage capacity | 22 |
|------------------|----|

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

### **Packaging**

|                        |               |
|------------------------|---------------|
| Minimum order qty.     | 1             |
| Piece(s) / inner box   | 1             |
| Piece (s) / outer case | 1             |
| Piece(s) / pallet      | 40            |
| Case volume            | 0.0052 cbm    |
| Weight of product      | 1.3595 kg     |
| Weight of case         | 1.36 kg       |
| EAN of inner box       | 4038917156257 |
| EAN case               | 4038917156240 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

