



Product description

| | |
|---------------------|--|
| Order number | 67.759.001 |
| Product description | UV cuvette, max. working volume: 2 ml, material: special plastic, transparent, (HxW): 45 x 12.5 mm, optical sides: 2, suitable for use beginning at 220 nm, path length: 10 mm, centre height: 15 mm, sorted by mould cavity number, free from DNA/RNase/protein, 100 piece(s)/polystyrene container |

Product characteristics

| | |
|--------------------------------|--------------------------------------|
| Closure type | push cap |
| Application | suitable for use beginning at 220 nm |
| Volume of work | 50 µl - 2 ml |
| Number of optical sides | 2 |
| Path length | 10 mm |
| Centre height | 15 mm |
| Wavelength in the cuvette area | 220 nm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|---------|
| Max. volume | 2.5 ml |
| Width of product | 12.5 mm |
| Height of product | 45 mm |

Material & colours

| | |
|-------------------|----------------------------------|
| Product material | special plastic |
| Colour of product | transparent |
| Closure material | Low Density Polyethylene (LD-PE) |
| Colour of cap | red |

Purity & certification

| | |
|------------------|-----------------------------|
| Product category | no medical device IVD |
| Purity standard | free from DNA/RNase/protein |
| Batched | yes |

Packaging

| | |
|-------------------------------|-----------------------|
| Minimum order qty. | 400 |
| Type of smallest subpackaging | polystyrene container |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 400 |
| Piece(s) / pallet | 31200 |
| Depth of box | 222 mm |
| Width of box | 202 mm |
| Height of box | 56 mm |
| Depth of case | 418 mm |
| Width of case | 278 mm |
| Height of case | 130 mm |
| Case volume | 0.0151 cbm |
| Weight of product | 0.0018 kg |
| Weight of case | 1.41 kg |
| EAN of inner box | 4038917264501 |
| EAN case | 4038917264495 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

