



Product description

Order number	95.987.002
Product description	PCR work tray, material: PC, transparent, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, 5 piece(s)/bag

Product characteristics

Suitable for	0.1 ml and 0.2 ml PCR tubes
Format	12 x 8

Size

Width of product	79 mm
Length of product	115 mm

Material & colours

Product material	Polycarbonate (PC)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
PCR work tray, PC, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 5
Type of smallest subpackaging bag
Piece(s) / inner box 5
Piece (s) / outer case 100
Piece(s) / pallet 32400
Depth of case 268 mm
Width of case 191 mm
Height of case 100 mm
Case volume 0.0051 cbm
Weight of product 0.0106 kg
Weight of case 1.25 kg
EAN of inner box 4038917124157
EAN case 4038917019118

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

