



Product description

Order number	95.988.003
Product description	Base station, material: PC, green, (LxWxH): 128 x 85.5 x 30.4 mm, suitable for work tray 95.987.002 & 0.2 ml PCR tubes, 5 piece(s)/case

Product characteristics

Suitable for	work tray 95.987.002 & 0.2 ml PCR tubes
--------------	---

Size

Width of product	85.5 mm
Height of product	30.4 mm
Length of product	128 mm

Material & colours

Product material	Polycarbonate (PC)
Colour of product	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	5
Type of smallest subpackaging	case
Piece(s) / inner box	5
Piece (s) / outer case	75
Piece(s) / pallet	3600
Depth of case	598 mm
Width of case	398 mm
Height of case	145 mm
Case volume	0.0345 cbm
Weight of product	0.0302 kg
Weight of case	2.83 kg
EAN of inner box	4038917124188
EAN case	4038917019149

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

