



Product description

| | |
|---------------------|--|
| Order number | 95.988.004 |
| Product description | Base station, material: PC, yellow, (LxWxH): 128 x 85.5 x 30.4 mm, suitable for work tray 95.987.002 & 0.2 ml PCR tubes, 5 piece(s)/case |

Product characteristics

| | |
|--------------|---|
| Suitable for | work tray 95.987.002 & 0.2 ml PCR tubes |
|--------------|---|

Size

| | |
|-------------------|---------|
| Width of product | 85.5 mm |
| Height of product | 30.4 mm |
| Length of product | 128 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polycarbonate (PC) |
| Colour of product | yellow |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 5
Type of smallest case
subpackaging
Piece(s) / inner box 5
Piece (s) / outer case 75
Piece(s) / pallet 3600
Depth of case 598 mm
Width of case 398 mm
Height of case 145 mm
Case volume 0.0345 cbm
Weight of product 0.0302 kg
Weight of case 2.83 kg
EAN of inner box 4038917124195
EAN case 4038917019156

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

