



Product Specification
PCR single tube, 0.5 ml, PCR Performance
Tested, colour mix, PP, flat cap

Page 1



Product description

Order number	72.735.992
Product description	PCR single tube, 0.5 ml, PCR Performance Tested, red/green/yellow/blue/violet/orange, material: PP, with anti-contamination shield, flat cap, per colour 500 piece(s)/minigrip bag

Product characteristics

Graduation	yes
Cap	flat cap
Closure type	push cap
Suitable for qPCR	no
Base shape	conical
High or low profile	High Profile
Number of wells	1
Volume of work	25 µl - 0,5 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
PCR single tube, 0.5 ml, PCR Performance
Tested, colour mix, PP, flat cap

Page 2

Size

Sample volume	0.5 ml
Max. volume	500 µl
Diameter	7.9 mm
Width of product	10 mm
Height of product	29.5 mm
Length of product	30.9 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	colour mix

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	3000
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	3000
Piece (s) / outer case	3000
Piece(s) / pallet	198000
Depth of case	428 mm
Width of case	323 mm
Height of case	155 mm
Case volume	0.0214 cbm
Weight of product	0.0004 kg
Weight of case	1.67 kg
EAN of inner box	4038917120067
EAN case	4038917017053

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT