



### Product description

Order number	72.985.092
Product description	PCR strip of 8, up to 200 µl, High Profile, PCR Performance Tested, white, material: PP, without cap, optimised for qPCR, 120 piece(s)/minigrip bag

### Product characteristics

Graduation	no
Closure type	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	High Profile
Number of wells	8
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Sample volume	0.2 ml
Max. volume	200 µl
Diameter	5.4 mm
Width of product	8.1 mm
Height of product	20.7 mm
Length of product	70.9 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	white

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

### **Packaging**

Minimum order qty.	480
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	120
Piece (s) / outer case	480
Piece(s) / pallet	155520
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0011 kg
Weight of case	0.67 kg
EAN of inner box	4038917374026
EAN case	4038917374019

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

