



**Product Specification**  
**PCR strip of 8, PCR Performance Tested,**  
**transparent, PP, round cap**

**Page 1**



**Product description**

Order number	72.994.002
Product description	PCR strip of 8, High Profile, PCR Performance Tested, transparent, material: PP, round cap, cap attached

**Product characteristics**

Graduation	no
Cap	round cap
Closure type	push cap
Base shape	conical
High or low profile	High Profile
Number of wells	8

**Size**

Sample volume	0.2 ml
---------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**PCR strip of 8, PCR Performance Tested,**  
**transparent, PP, round cap**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	240
Piece(s) / inner box	240
Piece (s) / outer case	1200
Piece(s) / pallet	1200
Case volume	0.0299 cbm
Weight of product	0.0031 kg
Weight of case	4.3 kg
EAN of inner box	4038917124324
EAN case	4038917080736

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

