




Product Specification
PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, Low DNA-binding, PCR Performance Tested, PP

Page 1



Product description

| | |
|---------------------|---|
| Order number | 72.1979.700 |
| Product description | PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, Low DNA-binding, PCR Performance Tested, material: PP, 10 piece(s)/minigrip bag |

Product characteristics

| | |
|---------------------------------|----------------|
| Suitable for qPCR | yes |
| Base shape | conical |
| Construction of the plate edges | half skirt |
| High or low profile | High Profile |
| Number of wells | 96 |
| Volume of work | 10 µl - 200 µl |
| Centrifugation max (RCF) | 3000 x g |
| Centrifugation temperature | 22 |
| Centrifugation time | 5 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
PCR plate half skirt, 96 well, transparent, High
Profile, 200 µl, Low DNA-binding, PCR
Performance Tested, PP

Page 2

Size

| | |
|-------------------|----------|
| Sample volume | 0.2 ml |
| Max. volume | 200 µl |
| Diameter | 4.95 mm |
| Width of product | 86 mm |
| Height of product | 19.33 mm |
| Length of product | 126.1 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 100 |
| Type of smallest subpackaging | minigrip bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 100 |
| Piece(s) / pallet | 5600 |
| Depth of box | 270 mm |
| Width of box | 230 mm |
| Height of box | 131 mm |
| Depth of case | 282 mm |
| Width of case | 280 mm |
| Height of case | 247 mm |
| Case volume | 0.0195 cbm |
| Weight of product | 0.0167 kg |
| Weight of case | 2.29 kg |
| EAN of inner box | 4038917574631 |
| EAN case | 4038917574624 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT