



**Product Specification**  
**PCR plate full skirt, 96 well, transparent, Low Profile, 100 µl, Biosphere® plus, PP**

**Page 1**



**Product description**

Order number	72.1980.201
Product description	PCR plate full skirt, 96 well, transparent, Low Profile, 100 µl, Biosphere® plus, material: PP, 1 piece(s)/blister

**Product characteristics**

Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	full skirt
High or low profile	Low Profile
Number of wells	96
Volume of work	5 µl - 100 µl
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**PCR plate full skirt, 96 well, transparent, Low**  
**Profile, 100 µl, Biosphere® plus, PP**

**Page 2**

**Size**

Sample volume	0.1 ml
Max. volume	100 µl
Diameter	5.5 mm
Width of product	85.5 mm
Height of product	15.5 mm
Length of product	127.8 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	20
Type of smallest subpackaging	blister
Piece(s) / inner box	10
Piece (s) / outer case	20
Piece(s) / pallet	1120
Depth of box	270 mm
Width of box	230 mm
Height of box	131 mm
Depth of case	282 mm
Width of case	280 mm
Height of case	247 mm
Case volume	0.0195 cbm
Weight of product	0.0174 kg
Weight of case	0.98 kg
EAN of inner box	4038917573238
EAN case	4038917573221

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**