




Product Specification
PCR plate full skirt, 96 well, transparent, Low Profile, 100 µl, Low DNA-binding, PCR Performance Tested, PP

Page 1



Product description

Order number	72.1980.700
Product description	PCR plate full skirt, 96 well, transparent, Low Profile, 100 µl, Low DNA-binding, PCR Performance Tested, material: PP, 10 piece(s)/minigrip bag

Product characteristics

Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	full skirt
High or low profile	Low Profile
Number of wells	96
Volume of work	5 µl - 100 µl
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
PCR plate full skirt, 96 well, transparent, Low
Profile, 100 µl, Low DNA-binding, PCR
Performance Tested, PP

Page 2

Size

Sample volume	0.1 ml
Max. volume	100 µl
Diameter	5.5 mm
Width of product	85.5 mm
Height of product	15.5 mm
Length of product	127.8 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	50
Piece (s) / outer case	100
Piece(s) / pallet	5600
Depth of box	270 mm
Width of box	230 mm
Height of box	131 mm
Depth of case	282 mm
Width of case	280 mm
Height of case	247 mm
Case volume	0.0195 cbm
Weight of product	0.0172 kg
Weight of case	2.34 kg
EAN of inner box	4038917574679
EAN case	4038917574662

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-24

This document was prepared by EDP support and is valid without signature.



SARSTEDT