



## Product description

|                     |   |
|---------------------|---|
| Order number        | 95.1999   |
| Product description | qPCR film, free from DNase/RNase, material: PO, highly transparent, suitable for qPCR, encapsulated adhesive, application area: -70 - 100 °C, heat sensitive adhesive, lowest evaporation rates, 100 piece(s)/bag |

## Product characteristics

|                     |        |
|---------------------|--------|
| Closure type        | film   |
| Light protection    | no     |
| Suitable for qPCR   | yes    |
| Minimum temperature | -70 °C |
| Maximum temperature | 100 °C |

## Size

|                   |        |
|-------------------|--------|
| Width of product  | 77 mm  |
| Length of product | 140 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polyolefin (PO)    |
| Colour of product | highly transparent |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Purity standard  | free from DNase/RNase   |
| Batched          | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 100           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1600          |
| Piece(s) / pallet             | 192000        |
| Depth of case                 | 173 mm        |
| Width of case                 | 135 mm        |
| Height of case                | 283 mm        |
| Case volume                   | 0.0066 cbm    |
| Weight of product             | 0.0021 kg     |
| Weight of case                | 3.92 kg       |
| EAN of inner box              | 4038917303590 |
| EAN case                      | 4038917388238 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

