



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 93.844.100   |
| Product description | Rack, material: PC, transparent, format: 10 x 2, (LxWxH): 257 x 62 x 55 mm, for 20 collection tubes, suitable for tubes, all S-Monovette® diameters, 1 piece(s)/case |

**Product characteristics**

|                        |                                   |
|------------------------|-----------------------------------|
| Suitable for           | tubes, all S-Monovette® diameters |
| Storage capacity       | 20                                |
| Format                 | 10 x 2                            |
| Autoclavable (max. °C) | 121 °C                            |

**Size**

|                   |        |
|-------------------|--------|
| Width of product  | 62 mm  |
| Height of product | 55 mm  |
| Length of product | 257 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polycarbonate (PC) |
| Colour of product | transparent        |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | no                      |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1             |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 42            |
| Piece(s) / pallet             | 1260          |
| Depth of box                  | 65 mm         |
| Width of box                  | 58 mm         |
| Height of box                 | 260 mm        |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 340 mm        |
| Case volume                   | 0.0539 cbm    |
| Weight of product             | 0.1907 kg     |
| Weight of case                | 10.14 kg      |
| EAN of inner box              | 4038917121712 |
| EAN case                      | 4038917017961 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

