



**Product Specification**  
**Block Rack D13, Ø opening: 13 mm, 5 x 10,  
orange**

**Page 1**



**Product description**

Order number	93.853.135
Product description	Block Rack D13, for 50 collection tubes, Ø opening: 13 mm, format: 5 x 10, orange, material: PP

**Product characteristics**

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

**Size**

Diameter of opening	13 mm
Width of product	96 mm
Height of product	45 mm
Length of product	180 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Block Rack D13, Ø opening: 13 mm, 5 x 10,  
orange**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	orange

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	16
Piece(s) / pallet	1440
Depth of case	398 mm
Width of case	368 mm
Height of case	120 mm
Case volume	0.0176 cbm
Weight of product	0.1687 kg
Weight of case	3.16 kg
EAN of inner box	4038917122566
EAN case	4038917074063

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**